DAO (Diamine Oxidase) REA

Call 800-592-5726 or [Request Information](#)

 For Research Use Only

Sample Type: serum, plasma

Sample Size: 100 µl

Size: 100 tests

Range: 2.1 - 80 U/ml

Sensitivity: 0.2 U/ml

Incubation Time: 2 hour(s) 30 min(s)

Protocol:  [DAO \(Diamine Oxidase\) REA](#)

International Distribution: Only available in North America

Price/Order:

[Buy Now](#)

[Price](#)

Availability: [Please Order In Advance](#)

Please call or email ALPCO
in advance for availability
800-592-5726

Related Products:

➤ [Histamine ELISA](#)
➤ [Hemoglobin / Haptoglobin ELISA](#)
➤ [S100A12 ELISA](#)

This item appears in: ➤ [Gastroenterology](#)

Summary:

Diamine oxidase (DAO) is a body's own enzyme that metabolizes histamine. Although DAO is found practically in the whole body, the most important site of its action is the intestine. The enzymatic activity of DAO determines the histamine degradation speed. In the case of DAO deficiency or inhibition, incorporated or endogenous histamine cannot be degraded quickly enough, and the symptoms of histamine intolerance are presented. Millions of people suffer from gastrointestinal problems, migraine, irritations of nasal mucosa and other allergy-like symptoms after consumption of certain nutrients. Too much histamine in the body can be the reason for this wide range of symptoms. Another possibility for reduced DAO function could be the intake of activity-inhibiting substances, such as alcohol or medication. Histamine induced food intolerance is not IgE-mediated. Determination of the DAO activity in serum or plasma is a suitable marker for diagnosis of histamine intolerance and the associated symptoms.

The described Radioextractionassay (REA) is intended for the quantitative determination of diamine oxidase activity in serum and plasma.

[Copyright © 2009 ALPCO Diagnostics](#) | [26-G Keewaydin Drive](#) | [Salem, NH 03079](#) | [Toll Free: 800-592-5726](#) | [Fax: 603-812-5726](#)
[Terms of Sale](#) | [Site Map](#) | [Feedback](#) | [Employee Login](#) | [Distributor Login](#)

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
Additional information about our products is available on our homepage

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15. LITERATUR

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). Agents and Actions 23:361-65.
2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. Scand J Clin Lab Invest 24:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. Inflammation Research 48:169-70.
4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. Allergy Proceedings 15:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. Inflammation Research 47:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millennium. Curr Probl Dermatol, Basel, Karger, 28:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches Clin Exp Allergy 23: 982-85.
8. Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel Allergologie 9: 426-30.
9. Jarisch R (1999) Histamin Intoleranz. Stuttgart, Thieme

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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Manual


DAO ELISA

For the in vitro determination of DAO in serum and plasma

Valid from 28.10.2010

K 8500

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of diamine oxidase (DAO) in serum and plasma.

2. INTRODUCTION

Diamine oxidase (DAO) is a body's own enzyme that metabolizes histamine. Although DAO is found practically in the whole body, the most important site of its action is the intestine. The enzymatic activity of DAO determines the histamine degradation speed. In the case of DAO deficiency or inhibition, incorporated or endogenous histamine cannot be degraded quickly enough, and the symptoms of histamine intolerance are presented. Millions of people suffer from gastrointestinal problems, migraine, irritations of nasal mucosa and other allergy-like symptoms after consumption of certain nutrients. Too much histamine in the body can be the reason for this wide range of symptoms.

Determination of the DAO concentration in serum and plasma is a suitable marker for diagnosis of histamine intolerance and associated symptoms.

Our DAO-ELISA kit is intended for determination of the diamine oxidase (DAO) concentration in serum and plasma.

Indications

- Frequent headaches or migraine
- Snuffles after consumption of histamine-containing nutrients
- Tissue oedema
- Eyelid turgor
- Skin redness
- Limb aches
- Gastrointestinal discomfort
- Monitoring of a histamine free diet

3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 8500MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 8500WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 8500ST	STD	Standards, lyophilized (see specification for concentration range)	4 x 7 vials
K 8500KO1	CTRL	Control, lyophilized (see specification for concentration range)	4 x 1 vial
K 8500KO2	CTRL	Control, lyophilized (see specification for concentration range)	4 x 1 vial
K 8500A2	AB	Detection antibody, (biotinylated), concentrate	1 x 200 µl
K 8500K	CONJ	Conjugate, (streptavidin peroxidase labeled), concentrate	1 x 200 µl
K 8500SV	STDBUF	Standard dilution buffer, ready to use	1 x 25 ml
K 8500PV	SAMPLEBUF	Sample dilution buffer, ready to use	1 x 50 ml
K 8500TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 8500AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Standard laboratory reaction vessels (1.5 ml)
- Standard laboratory reaction vessel (15 ml)
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Vortex mixer
- Bidistilled water (aqua bidest.)
- Shaking incubator at 37°C
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. The **WASHBUF** (wash buffer concentrate) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The lyophilized **STD** (standards) and **CTRL** (controls) are stable at 2-8°C until the expiry date stated on the label. The **STD** (standards) and **CTRL** (controls) must be reconstituted with **500 µl STDBUF** (standard dilution buffer). Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards are **not stable**.
- The **AB** (detection antibody, biotinylated) must be diluted **1:101** in wash buffer (e.g. 100 µl AB + 10 ml wash buffer). The undiluted AB is stable at **2-8 °C** until the expiry date given on the label. Diluted antibody solution is **not stable** and can not be stored.
- The **CONJ** (conjugate, POD-antibody) must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted CONJ is stable at **2-8 °C** until the expiry date stated on the label. Diluted conjugate is **not stable** and can not be stored.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

6. SAMPLE PREPARATION

Serum / Plasma

Pipet **25 µl** of fresh sample in a 1,5 ml reaction vial, add **225 µl SAMPLEBUF** (sample dilution buffer) and mix well (corresponds to **1:10 dilution**).

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the “sandwich” technique with two polyclonal antibodies against recombinant DAO.

Standards, control and diluted samples which are assayed for DAO are added into the wells of a micro plate coated with polyclonal rabbit anti-DAO antibody. During the first incubation step, DAO is bound by the immobilized primary antibody. Then a biotinylated polyclonal anti-DAO antibody, is added into each microtiter well. In the next step, the streptavidin-POD-conjugate is added and a “sandwich” of

1st antibody – DAO - biotinylated antibody – streptavidin-POD-conjugate

is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of DAO. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. standard concentration is generated, using the values obtained from the standard.

Test procedure

1.	Prior to use in the assay allow all reagents and samples to come to room temperature and mix by gentle swirling and inversion.
2.	Mark the positions of STD (Standards)/ SAMPLE/CTRL (Controls) in duplicate on a protocol sheet
3.	Take PLATE (microtiter strips) out of the kit. Store unused strips in the original package bag at 2-8° C. Strips are stable until expiry date stated on the label
4.	Wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
5.	Add 100 µl of STD (Standards)/ SAMPLE/CTRL (Controls) in duplicate into respective well
6.	Cover the plate tightly and incubate for 2 hours at 37° C on a horizontal mixer
7.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
8.	Add 100 µl of AB (Detection antibody, 2nd Antibody) into each wells mix gently
9.	Cover the plate tightly and incubate for 1 hour at 37° C on a horizontal mixer

10.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
11.	Add 100 µl of CONJ (Conjugate) into each well
12.	Cover the plate tightly and incubate for 1 hour at 37° C on a horizontal mixer
13.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
14.	Add 100 µl of SUB (Substrate) into each well
15.	Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark
16.	Add 50 µl of STOP (Stop solution) into each well, shake well
17.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

*The intensity of the colour change is temperature sensitive. We recommend to observe the colour change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-Parameter-algorithm".

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Samples

The obtained DAO concentration must be multiplied by a factor of **10**.

9. LIMITATIONS

Samples with DAO concentrations **outside the standard curve** range should be diluted further with **SAMPLEBUF** (sample dilution buffer) to obtain readings within the standard curve and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal range

< 3 U/ml:	high incidence for HIT (Histamine intolerance)
3 - 10 U/ml:	HIT probable
> 10 U/ml:	low incidence for HIT

Normal concentration range: Jarisch R et al. (1999)

It is recommended for each laboratory to establish its own normal range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) was calculated from 16 replicate determinations on each of two samples.

Intra-Assay CV n= 16

Sample	DAO [U/ml]	Intra-Assay CV [%]
1	5.0	1.42
2	5.9	1.72

The total precision (inter-assay variation) was calculated from data on 2 samples obtained in 8 different assays by different technicians on three different days.

Inter-Assay CV n= 8

Sample	DAO [U/ml]	Inter-Assay CV [%]
1	23.27	7.9
2	13.36	10.7

Sensitivity

The sensitivity limit was set as $B_0 + 3 \text{ SD}$. The Zero-standard was measured 20 times.

Sensitivity n=20

Sample	Mean value [OD]	Standard variation	Detection limit [U/ml]
1	0.019	0.0036	0.52

Sample dilution

Two patient samples were diluted and assayed. The results are shown below:

Linearity n= 2

Sample	Dilution	Expected [U/ml]	Measured [U/ml]
A	undiluted	15,65	15,65
	1:10	7,83	7,68
	1:20	3,92	3,86
	1:40	1,96	1,86
B	undiluted	5,58	5,58
	1:10	2,79	2,89
	1:20	1,39	1,45
	1:40	0,7	0,7

12. PRECAUTIONS

- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulphuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

15. REFERENCES

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). Agents and Actions 23:361-65.
2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. Scand J Clin Lab Invest 24:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. Inflammation Research 48:169-70.
4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. Allergy Proceedings 15:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. Inflammation Research 47:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millenium. Curr Probl Dermatol, Basel, Karger, 28:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches Clin Exp Allergy 23: 982-85.
8. Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel Allergologie 9: 426-30.
9. Jarisch R (1999) Histamin Intoleranz. Stuttgart, Thieme

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

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
Radioextractionassay for the quantitative determination of
diamine oxidase activity in serum

Gültig ab / Valid from 09.06.2008

K 8220

Immundiagnostik AG

For Research Use Only
For Reference Purposes Only

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13. LITERATUR

- Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). Agents and Actions 23:361-65.
- Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. Scand J Clin Lab Invest 24:163-68.
- Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. Inflammation Research 48:169-70.
- Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model for food intolerance. Allergy Proceedings 15:27-32.
- Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. Inflammation Research 47:396-400.
- Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millennium. Curr Probl Dermatol, Basel, Karger, 28:64-73.
- Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches Clin Exp Allergy 23: 982-85.

- Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel Allergologie 9: 426-30.
- Jarisch R (1999) Histamin Intoleranz. Stuttgart, Thieme

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
DAO-REA (^3H)

Radioextraction assay for the quantitative determination
of diamine oxidase activity in serum

Valid from 09.06.2008

K 8220

Immundiagnostik AG

Content

1. INTENDED USE	13
2. INTRODUCTION	13
3. MATERIAL SUPPLIED	14
4. MATERIAL REQUIRED BUT NOT SUPPLIED	14
5. STORAGE OF REAGENTS	15
6. ASSAY PROCEDURE	15
Principle of the test	15
Procedural notes	15
Preparation	15
Test procedure	16
7. RESULTS	17
8. QUALITY CONTROL	17
Expected values	17
9. PERFORMANCE CHARACTERISTICS	18
Precision and reproducibility	18
Sensitivity	18
Recovery:	19
10. PRECAUTIONS	19
11. TECHNICAL HINTS	19
12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	20
13. REFERENCES	20

1. INTENDED USE

The described Radioextractionassay (REA) is intended for the quantitative determination of diamine oxidase activity in serum and plasma. It is for in vitro diagnostic use only.

2. INTRODUCTION

Diamine oxidase (DAO) is a body's own enzyme that metabolizes histamine. Although DAO is found practically in the whole body, the most important site of its action is the intestine. The enzymatic activity of DAO determines the histamine degradation speed. In the case of DAO deficiency or inhibition, incorporated or endogenous histamine cannot be degraded quickly enough, and the symptoms of histamine intolerance are presented. Millions of people suffer from gastrointestinal problems, migraine, irritations of nasal mucosa and other allergy-like symptoms after consumption of certain nutrients. Too much histamine in the body can be the reason for this wide range of symptoms. Another possibility for reduced DAO function could be the intake of activity-inhibiting substances, such as alcohol or medication.

Histamine induced food intolerance is not IgE-mediated. Determination of the DAO activity in serum or plasma is a suitable marker for diagnosis of histamine intolerance and the associated symptoms.

With our easy-to-use, reliable and standardised test kit it is possible to quantify the biological activity of DAO in the circulation. Only 100µl of serum is needed for the test, results are available within 3 hours.

Indications:

- Detection of histamine intolerance
- Monitoring of a histamine-free diet

3. MATERIAL SUPPLIED

Catalogue No.	Content	Kit Components	Quantity
K8220ST	STD	Standards with diamine oxidase from porcine kidney, lyophilized.	1 ml
K8220KO	CTRL	Diamine oxidase from porcine kidney, lyophilized	1 ml
K8220AP	ASYBUF	assay buffer, ready to use	15 ml
K8220SUB	SUB	radio-labelled putrescine-dihydrochloride in a stabilised solution., ready to use	5,5 ml
K8220EXT	EXTSOL	non-toxic, chlorine free extraction solution, ready to use	105 ml
K8220SZI	SZIN	Aquasafe Liquid Szintillation cocktail, ready to use	105 ml

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5. STORAGE OF REAGENTS

- Freshly collected samples may be stored at room temperature or at 2-8°C up to 24 hours. For prolonged storage samples must be frozen at -20°C.
- Avoid more than 3 freeze cycles.
- The reconstituted standards (STD) and the control (CTRL) are stable for 4 weeks at 2-8°C.
- All reagents are stable at 2-8°C until the expiry date stated on the label

6. ASSAY PROCEDURE

Principle of the test

Human serum may be used as an analyte. The activity of the diamine oxidase is determined by quantitating the reaction product. Radiolabelled putrescine-dihydrochloride is used as a substrate. The resulting Δ^1 pyrroline, containing the radiolabel, is extracted selectively from the matrix by a liquid extraction step. A non-toxic, chlorine free solvent with high capacity is used for the extraction.

Finally scintillation fluid is added to the organic phase containing the radiolabelled Δ^1 pyrroline and radioactivity is determined in a beta-counter. The signal is directly proportional to the activity of DAO in the sample.

Procedural notes

- Samples should be mixed well before assaying.
- We recommend duplicates for all values.

Preparation

- The lyophilized standards (STD) and the control (CTRL) must each be reconstituted with 1ml distilled H_2O . Allow the vial content to dissolve for 10 minutes; afterwards vortex well.
- All reagents should have reached room temperature (18°C to 26°C) before use. Prepare and label PP-tubes for standards, control and samples as appropriate.

Test procedure

Pipette 100 μl of reconstituted standards (STD), control (CTRL) and samples into the respective tubes
Add 100 μl of assay buffer (ASYBUF) into all tubes
Add 50 μl of substrate (SUB) to all tubes
Vortex all tubes thoroughly
Incubate all tubes for 150 min at 37°C , shake if possible
Add 1 ml of extraction solution (EXTSOL) into all tubes, mix thoroughly on a vortex mixer for at least 15 seconds This step is essential for a complete extraction!
Centrifuge the extraction mixture for 1 min at $> 3000 \times g$
Transfer 800 μl of the upper phase (organic) into a scintillation vial (Alternative method of running the test, see below *)
Add 1 ml of szintillaton (SZIN) fluid into all szinti-vials, mix well
Count every vial for 1 min in the beta counter

*Alternative method of running the test:

Decantation of the supernatant:

Alternatively the reaction tubes can be incubated at -20°C for at least 90 minutes to freeze the lower (red) phase. The organic supernatant remains liquid and can easily be decanted directly into the scintillation vials. It's necessary that the entire lower phase is frozen to make sure that it completely remains in the tube while decanting the organic phase.

Please avoid the usage of polystyrene tube racks or increase the incubation time at -20°C . Metal or plastic tube racks are recommended. A preincubation of these racks at -20°C can reduce the required incubation time. In any case it is highly recommended to make sure that a constant phase separation is warranted.

Transfer 800 μl of the upper phase (organic) into a scintillation vial
Add 1 ml of szintillaton (SZIN) fluid into all szinti-vials, mix well
Count every vial for 1 min in the beta counter

7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the „4-Parameter-algorithm“.

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e.g. 0.01).

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3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e.g. 0.01).

8. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

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Norm concentration range: Jarisch R et al. (1999)

We recommend each laboratory to establish its own norm concentration range.

A diagnosis of histamine-intolerance must be correlated to the actual concentration of histamine.
It is IMPOSSIBLE to state histamine intolerance in case of anaphylactic shock and during pregnancy!

9. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n=10)		
Sample	DAO-Activity (U/ml)	VK [%]
1	41,5	4,07
2	2,8	6,29
3	27,5	4,08

Inter-Assay (n=5)		
Sample	DAO-Activity (U/ml)	VK [%]
1	42,4	5,59
2	2,7	16,16
3	27,4	4,39

Sensitivity

The sensitivity was set as 0,2 U/ml (B0 + 3 SD). The zero-standard was measured 10 times.

Recovery:

2 serums were spiked with 10 U/ml DAO and measured with the assay.

Serum	Unspiked	Expected	Measured	Recovery
1	19,1	29,1 U/ml	29,7 U / ml	102,1 %
2	1,5	11,5 U/ml	9,6 U / ml	83,5 %

10. PRECAUTIONS

- For in vitro diagnostic use only.
- The quality control guidelines should be followed.
- Avoid all contact with the reagents by using gloves.
- Radioactive waste must be disposed according to local regulations.
- The extraction solvent is a non-halogenated organic solvent and has to be disposed according to local regulations.

11. TECHNICAL HINTS

- Extraction time must not be shorter than 15 sec as extraction recovery may vary with shorter extraction times.
- Do not use reagents beyond expiration date.
- Do not expose reagents to direct sunlight.
- Do not mix or substitute reagents with those from other lots or sources.
- The assay should always be performed according to the enclosed manual.

12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be sent to Immundiagnostik AG together with a written complaint.

13. REFERENCES

- Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). Agents and Actions 23:361-65.
- Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. Scand J Clin Lab Invest 24:163-68.
- Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. Inflammation Research 48:169-70.
- Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. Allergy Proceedings 15:27-32.
- Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. Inflammation Research 47:396-400.
- Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millenium. Curr Probl Dermatol, Basel, Karger, 28:64-73.
- Wantke F et al. (1993) Histamine free diet: treatment of choice for histamine induced food intolerance and supporting treatment for chronic headaches Clin Exp Allergy 23: 982-85.
- Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel Allergologie 9: 426-30.
- Jarisch R (1999) Histamin Intoleranz. Stuttgart, Thieme

Used Symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Radioactivity

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